

## **Dynamic weight bearing 動態負重測試**

### **1. Purpose**

- 1.1 To examine the spontaneous pain of a mouse by using bioseb dynamic weight bearing system.

### **2. Safety Requirements**

- 2.1 General laboratory procedures should be followed, which include: no eating, no chewing gum, no drinking, and no applying of cosmetics in the work area. The researcher must wear the laboratory coat, gloves and a mask during the experiment.

### **3. Associated Documents**

- 3.1 [https://www.bioseb.com/bioseb/anglais/default/item\\_id=1213\\_cat\\_id=5\\_Advanced+Dynamic+Weight+Bearing.php](https://www.bioseb.com/bioseb/anglais/default/item_id=1213_cat_id=5_Advanced+Dynamic+Weight+Bearing.php)
- 3.2 <https://www.sunpointworld.com/products/dynamic-weight-bearing/>

### **4. Notes**

- 4.1 The mice must be maintained in a controlled environment with stable temperature, humidity, and air pressure. Keep minimal disturbing and stress on the tested mice whenever is possible.
- 4.2 The majority of mouse behavioral studies are age/sex/strain dependent. It is important to have the age-matched mice with the same sex and strain background in a single experiment.
- 4.3 The dynamic weight bearing measures the balance of mice limbs and tail. Therefore, mice toes should not be injured and tail cut (for genotyping) should be minimal.

### **5. Quality Control**

- 5.1 The mouse shall be habituated in the testing room 30 min before the test.
- 5.2 Calibrate the image and pressure sensor of dynamic weight bearing system before the starting of the measurement.
- 5.3 Record the time whenever the system shows detected error.

### **6. Equipment**

- 6.1 A computer for dynamic weight bearing.
- 6.2 An image capture of dynamic weight bearing.
- 6.3 A pressure sensor of dynamic weight bearing.
- 6.4 A cage of dynamic weight bearing.
- 6.5 A partition to reduce the environmental noise.

### **7. Supplies**

- 7.1 Cap.
- 7.2 Gloves.
- 7.3 Facial Mask.
- 7.4 Ethanol 70%.
- 7.5 Paper towel.
- 7.6 Kimwipes.
- 7.7 HOCl.

## **8. Procedures**

- 8.1 To measure the body weight of each mouse.
- 8.2 The mouse shall be habituated in the room for 30 min before the test.
- 8.3 Calibrate the image and pressure sensor of dynamic weight bearing system.
- 8.4 Enter the information of body weight and set up the system.
- 8.5 Put the mouse in the cage of dynamic weight bearing system and record for 5 minutes.
- 8.6 Place the mouse back to the home cage.
- 8.7 Repeat steps 1-4 until all mice finish the dynamic weight bearing test.
- 8.8 Analysis the data and export the data.

